#1 CHOICE OF HOSPITALS AND MOTHERS

Harmony™

MANUAL BREASTPUMP

BOMBA DE LACTANCIA MANUAL

TIRE-LAIT MANUEL

**Instructions
& Assembly**

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**Instructions
et Assemblage**


Table of Contents-English

Harmony Parts & Assembly	2-3
Ordering Parts	2
Before You Begin	4
Prior to Each Use	4
Hygiene	4
Single User Caution	4
Harmony Disassembly	5-6
Harmony Assembly	6-7
Using Your Breastpump	8
When to Use	9
Upgrading to electric pumping	9-12
Tips for Success & Trouble Shooting	12-13
Storing & Handling Breastmilk	13-14

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
The Harmony Breastpump.... from Medela®

The Harmony is a breakthrough in manual breastpump technology. It is the first and only manual pump with 2-Phase Expression® – a dual pumping pattern that comes closer than ever to the way a baby naturally breastfeeds.

The first pattern, called the Stimulation phase, mimics a baby's initial rapid, shallow sucking that encourages fast let-down. The second pattern, the Expression phase, simulates a baby's more relaxed, deeper sucking that promotes efficient milk removal.

Your Harmony Breastpump is designed to give you gentle comfort. The contoured handle is ergonomically designed to fit either hand comfortably. The handle also swivels, allowing you to position your hand and the pump more easily and comfortably on each breast. The SoftFit™ Breastshield is designed to gently massage the breast while you pump. As a one-piece unit, there is no assembly or extra pieces to clean and store.

Harmony....blending quiet,
effective pumping with ease,
gentle comfort and convenience.

Harmony Breastpump Parts & Assembly


A. Handle.....	#8007170
A1. Let-Down Symbol	
B. Diaphragm.....	#8207043
B1. Stem w/ O-Ring.....	Stem #8107119
	O-Ring #8247018
C. SoftFit Breastshield.....	#8100402
D. Pump Body.....	#8007196
E. White Membrane.....	#8207013
F. Yellow Valve.....	#8100475
G. Breastmilk Container.....	#8117009
H. Container Stand.....	#8100462
J1. Solid Container Lid.....	#8107057
J2. Container Lid with Hole.....	#8100473
J3. Disc Insert.....	#8107062

Note: Handle (A) & Diaphragm/Stem (B) do not need to be disassembled and cleaned after each use, unless it is apparent the assembly has become soiled. White O-Ring on stem does not need to be removed.

Ordering Parts

In the US and Canada, if your retailer does not carry the parts you need, you can place a credit card order by phone by calling 1-800-435-8316 (US & Canada), 1-815-363-1166 (US only) or 1-905-608-7272 (Canada). VISA, MasterCard, American Express and Discover Card are accepted. In other countries, please contact your local distributor. Please visit our website, www.medela.com for a list of worldwide distributors.

Breastpump Parts

Before You Begin

Before using your Harmony for the first time:

The Harmony sanitized version does not require washing before first use.

Disassemble parts – refer to diagram on page 3 and photos on page 5-6. Sanitize by boiling for 10 minutes.

Prior to each use:

Disassemble the SoftFit breastshield, pump body, valve, membrane and container.

Under normal use, the Pump Handle (A) and Diaphragm/Stem (B) only need to be cleaned periodically or when it is apparent they have become soiled. Hint: Do not remove small white O-ring from stem.

Wash remaining parts in dishwasher or in soapy water and rinse well. Air dry parts on clean towel and cover when not in use.

Do not use with Quick Clean™ Micro-Steam™ bags. Wash only according to the instructions above.

Hygiene:

Wash hands thoroughly before pumping. If your physician recommends it, clean breast and surrounding area with a damp cloth. Do not use soap or alcohol, which can dry the skin and cause cracked nipples.

Using a breastpump should never be uncomfortable. If you experience any pain, contact your healthcare provider or lactation specialist immediately.

Single User Caution

The Harmony Breastpump is a personal care item that should not be shared between users.

Harmony Disassembly

1. Remove the SoftFit breast shield (C) from the pump body (D).

2. Remove the pump handle assembly from the pump body (A-B). Check to see if it needs cleaning.

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
4. Pull yellow valve (F) from pump body (D) by twisting gently.



Harmony Disassembly (continued)

5. Remove white membrane (E) from yellow valve (F).

Harmony Assembly

1. Snap white membrane (E) button into small opening on yellow valve (F) until membrane lies completely flat on valve. (Both parts should be completely dry before assembly.)
2. Push yellow valve securely onto inner ring of pump body (D).
3. Screw collection container (G) onto pump body.

4. If the handle components (A & B) have been disassembled: place the oval shaped base of the stem (B1) into the diaphragm (B).

Hint: Make sure the oval shape is lined up correctly in the diaphragm.

Correct

Incorrect

5. Place silicone diaphragm/stem (B) into pump handle (A) and press down firmly around entire edge to ensure proper seal.
6. Insert stem of handle assembly into opening on top of pump body (D) and push until you feel it snap in.
7. Insert SoftFit breastshield (C) into top opening on pump body (D).

Using Your Harmony Breastpump

1. Center breastshield over nipple. Caution: pushing shield too firmly into breast may actually reduce milk flow.
2. The adjustable handle can be positioned for your maximum pumping comfort. The level of vacuum is controlled by the compression and the speed of the handle stroke. You decide the pumping rhythm that is best for you.
3. Start pumping in the Stimulation mode, by pressing down on the milk drops symbol (♻️) on top of the pump handle. Press and release quickly until you see milk flow begin. This mode mimics baby's initial rapid sucking and is designed to stimulate milk flow/let-down.
4. Once let-down has begun, switch to the Expression mode by pressing the bottom of the pump handle. This will give you a longer stroke that mimics babies' deeper, more relaxed sucking and is designed to maximize milk flow.

5. Once you are finished pumping, parts should be cleaned following directions on page 4.

When to Use the Harmony Breastpump

The Harmony Breastpump...

- allows you to pump manually on an occasional basis. That is, you are available to nurse your baby for most feedings each day, except for a few occasions during the week.
- is ideal to collect breastmilk for a feeding that will be missed due to an evening out.
- is also ideal to relieve occasional breast engorgement or sore nipples.


When long-term pumping is required, e.g. when the baby is not available to feed at the breast at least part of the time each day (a premature infant, for example), it is recommended that you use a full-size hospital electric breastpump such as the Symphony®, Lactina® or Classic™ Breastpump. Each of these pumps may be rented for short to extended periods of time, depending upon your needs. For the names and locations of Medela electric breastpump rentals, please call 1-800-TELL YOU (1-800-835-5968) in the US or Canada.

Upgrading to Electric Pumping with Medela Symphony or Lactina Breastpumps.

Symphony Breastpump Setup

Requires membrane cap #8100667, membrane #820094 and tubing #8007179 or #8007180. Sold separately.

1. Remove handle assembly from Harmony Breastpump.


2. Insert end of tubing into opening on top of pump.


3. Assemble protective membrane and membrane cap.

4. Insert yellow end of tubing into port on membrane cap.

5. Push membrane cap firmly into place on pumping mechanism. Adjust for comfort once pumping begins.

Lactina Breastpump Setup

Requires kit #67090 and tubing #8007194 or #8007180, sold separately. See Medela's instruction booklet for detailed kit assembly and cleaning.


1. Remove handle assembly from Harmony Breastpump.


2. Insert end of tubing into opening on top of Harmony pump.

3. Slide piston into cylinder.


4. Screw cylinder tightly into pump connector.

5. Attach pump connector to Lactina, and rotate until tab is secure.

6. Push piston into pumping arm, and rotate until snug in clamp.

7. Insert other end of tubing into port. Insert plug into unused Lactina port.

8. Set the Lactina vacuum regulator to “MIN.” Adjust for comfort once pumping begins.

Tips for Success and Troubleshooting:

1. Prior to pumping, relax, take a few deep breaths and think of your baby.
2. Using a breastpump should be easy and comfortable. The first few times may seem a bit awkward, but you will quickly find your best position and pumping speed.
3. Each mother reaches let-down at a different pace. Sometimes, you may not need to use the Stimulation mode at all and can start right in with normal Expression pumping.
4. If you feel the pump is not giving sufficient suction, you can check several things:
 - a. Make sure the white membrane (E) lies completely flat in place on the yellow valve (F) and is clean.
 - b. Make sure the stem assembly (B1) on the pump handle is seated all the way into the pump body (D) until it will go no farther. Make sure the white o-ring is in place in the stem handle (B1).
 - c. Make sure the pump body (D) is firmly screwed onto the container (G) and that the breastshield (C) is firmly seated into the pump body (D).
 - d. Make sure the oval base of the stem (B1) is properly seated in the oval base of the diaphragm (B).
 - e. Finally, make sure the diaphragm (B) is pushed all

the way onto the circular rim on the underside of the handle (A).

- f. If assistance is still needed after reviewing this booklet, please contact Medela Customer Service at 1-800-435-8316, or 1-815-363-1166, Monday-Friday, 7:30 AM to 8:00 PM Central Time.

Storing and Handling Breastmilk

Storage:

(Please check with your hospital for specific storage instructions, if needed.)

1. Wash hands before touching containers and avoid touching the interior of containers or caps.
2. Pump or express into clean collection containers.
3. If you do not intend to use the expressed breastmilk within a few days, freeze it.
4. Label containers with the collection date if you are not using the milk immediately.
5. Freeze milk in small portions, 2 to 4 ounces per container. Smaller amounts thaw more quickly.
6. When freezing milk, do not fill containers more than 3/4 full to allow room for liquid to expand. Cool milk prior to freezing. Never refreeze breastmilk that has been frozen.
7. Use solid lids and seal container tightly if you intend to refrigerate breastmilk for several days or freeze it.
8. Refer to the following chart for maximum storage periods under varying conditions.

Breastmilk Storage Guidelines

Breastmilk Storage (For Healthy Term Babies)					
	Room Temperature	Cooler with 3 Frozen Ice Packs	Refrigerator	Self-contained Refrigerator Freezer Unit	Deep Freezer
Freshly expressed breastmilk	4 hours at 66-72°F (19-22°C) ¹	24 hours at 59°F (15°C) ¹	5-7 days at 32-39°F (0°-4°C) ²	3-4 months ³	6-12 months at 0°F (-19°C) ³
Thawed breastmilk (previously frozen)	Do not store	Do not store	24 hours ⁴	Never refreeze thawed milk	Never refreeze thawed milk

¹ Hamosh M, Ellis L, Pollock D, Henderson T, and Hamosh P: Pediatrics, Vol. 97, No. 4, April 1996. pp 492-497. (4 hours at 77° F/25° C).

² Sosa, Roberto; Barness, Lewis: AJDC, Vol. 141, Jan. 1987.

³ Lawrence R, and Lawrence R: Breastfeeding: A Guide For the Medical Profession, 1999, p.894.

⁴ Arnold L; Recommendations for Collection, Storage and Handling of a Mother's Milk for Her Own Infant in the Hospital Setting, 3rd Edition. The Human Milk Banking Association of North America, Inc. 1999, p. 18.

Warming/Thawing Breastmilk

1. Place sealed container of frozen or refrigerated breastmilk in a bowl of warm water for about 30 minutes. Or hold under warm running water.
2. Do NOT thaw breastmilk in a microwave oven. High temperatures can destroy important nutrients in the breastmilk and there is a danger of scalding baby.
3. Once thawed or warmed, swirl container to blend fat, which may have separated during storage. Test for temperature on the inside of your wrist. It should be lukewarm.
4. Feed thawed milk immediately or store in refrigerator for up to 24 hours. Do not let thawed breastmilk sit out at room temperature.
5. If you are traveling or transporting breastmilk to the baby's care giver, a cooler carrier equipped with cooling gel packs, will help keep breastmilk cool until you reach your destination.

For a breastfeeding specialist or to locate Medela rental or retail locations, call the Breastfeeding National Network™ (BNN™) at 1-800-TELL YOU; or visit the Medela website: www.medela.com.

The Medela Customer Service Team can be reached at 1-800-435-8316.

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⁴ Arnold L: Recommendations for Collection, Storage and Handling of a Mother's Milk for Her Own Infant in the Hospital Setting, 3rd Edition. The Human Milk Banking Association of North America, Inc. 1999, p. 18.

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Fax: 1-905-608-8720

Email: info@medela.ca

Medela, Lactina, and Symphony are registered trademarks of Medela. Harmony, Natural Expression, SoftFit, QuickClean, PureLan, Classic, Breastfeeding National Network and BNN are trademarks of Medela.

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Medela, Inc., P.O.
Medela Canada, Inc.,