#1 CHOICE OF HOSPITALS AND MOTHERS

Swing™

BREASTPUMP

BOMBA DE LACTANCIA

TIRE-LAIT

www.medela.com

Medela, Inc.,
P.O. Box 660, McHenry, IL
60051-0660, USA
Ph/Tel: 1-800-435-8316
Fax: 1-815-363-1246
customer.service@medela.com

Medela Canada, Inc.,
4090B Sladeview Crescent
Unit 2, Mississauga, Ontario,
Canada L5L 5Y5
Ph/Tel: 1-800-435-8316
Fax: 1-800-995-7867
info@medela.ca

Printed in the U.S.A. Medela, 2-Phase Expression, Symphony and Pump In Style are registered trademarks of Medela, Inc. in the USA and a trademark in other countries. Swing, Harmony, Maximum Comfort Vacuum, SoftFit, PersonalFit, Quick Clean, Micro-Steam, DayCarrier, CSF, PureLan, Breastfeeding National Network, and BNN are trademarks of Medela, Inc. in the USA and in other countries. Application for registration of Medela is pending in other countries. NursingStool is a registered trademark of Ortho Products. Patents Pending. All rights reserved.

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Congratulations

on your new Swing breastpump!

Swing is small and lightweight, and easy to carry anywhere. Swing is a single electric pump ideal for:

- collecting breastmilk while away from your baby
- maintaining a milk supply if breastfeeding is temporarily interrupted
- relieving engorgement by single pumping with simultaneous massage

Swing™ and

When it comes to breastfeeding, babies know best.

So, how do babies nurse?

Babies nurse in two distinct phases:

- Fast and light sucking to start your milk flowing
- A slower, deeper suck to get as much milk as quickly as possible

Swing contains 2-Phase Expression®, Medela's latest research-based technology that closely mimics a baby's actual nursing:

- Phase 1, the stimulation phase, is a rapid rhythm to get your milk flowing
- Phase 2, the expression phase, is a slower rhythm for gentle and efficient expression

The combination of these two phases leads to faster let down and milk flow.

2-Phase Expression® is found only in these Medela® pumps:

Harmony™

Swing™

Symphony™

(available through rental)

Pump In Style® Advanced
The Metro Bag™

Table of Contents - English

Swing Parts Guide	1
Swing Motor Guide	1
Important Safeguards	2
Assembly Instructions	3-4
Battery Operation	4
Pumping Instructions	5
Maximum Comfort Vacuum™	6
Mobility	7
Cleaning & Hygiene	7
Prior to First Use	8
After Each use	8
Tubing Care	8-9
Storing Breastmilk	9-10
Troubleshooting	10-11
Warranty	12
Ordering Parts	13
Breastfeeding Accessories	13

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
Swing motor guide

Lights up = On
Does not light up = Off

Slow flashing = stimulation phase
Constantly on = expression phase

Rapid flashing = Indicates vacuum at minimum or maximum

Medical Equipment
With Respect to Electrical Shock, Fire and
Mechanical Hazards Only In Accordance With
IEC 60601 -1 CAN/CSA 22.2 No. 601.1 22PL

IMPORTANT SAFEGUARDS

This is a single user product. Use by more than one person may present a health risk and performance questions, and voids the warranty.

When using electrical products, especially when children are present, basic safety precautions should always be followed, including the following:

READ ALL INSTRUCTIONS

DANGER: To reduce the risk of electrocution:

1. Always unplug electrical products immediately after using.
2. Do not use while bathing.
3. Do not place or store product where it can fall or be pulled into a tub or sink.
4. Do not place in or drop into water or other liquid.
5. Do not reach for a product if it has fallen into water. Unplug immediately.

WARNING: To reduce the risk of burns, electrocution, fire or injury to persons:

1. A product should never be left unattended when plugged into an electric outlet.
2. Close supervision is necessary when this product is used by, on, or near children or invalids.
3. Use the product only for its intended use as described in these instructions. Do not use attachments not recommended by the manufacturer.
4. Never operate this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water.
5. Keep the cord away from heated surfaces.
6. Do not use outdoors or operate where aerosol (spray) products are being used or where oxygen is being administered.
7. Plug the power cord into the Swing breastpump first and then into wall socket.
8. Only use the power cord that comes with the Swing breastpump.
9. This product is a personal care item and should not be shared among mothers.
10. Make sure the voltage of the power cord is compatible with the power source.
11. Never use while sleeping or when drowsy.
12. Wash the appropriate pump kit parts before each use.
13. If the suction feels uncomfortable or painful, decrease the suction level. If this doesn't help, turn product off and remove from your breast. Consult your healthcare professional.

SAVE THESE INSTRUCTIONS

Assembly Instructions - Swing Breastpump

1. Make sure breastshield and connector are securely attached.

2. Snap membrane onto valve until membrane lies completely flat.

Be sure to take apart membrane from valve when cleaning, and dry completely before reassembling.

3. Push assembled valve & membrane onto connector.

4. Screw container into connector.

5. Insert the round, yellow end of tubing into opening on the pump body.

6. Insert the triangle, yellow end of tubing into the triangle port on the motor.

7. Plug the metal end of the power cord into the AC adapter port as shown. Plug the cord into a standard electrical outlet.

Battery Operation

Requires four (4) AA batteries (sold separately).

Fresh batteries will provide approximately 2 to 4 hours of pumping time.

1. To install batteries, push the button towards the clip to open the cover of the battery compartment.

2. Insert batteries as indicated inside the battery compartment, and snap the cover back into place.

Tip: Once you are finished pumping, unplug tubing from back of breastshield prior to setting down bottle.

Pumping Instructions

Be comfortably seated and relax. Using a footstool, like the Medela NursingStool®, helps you maintain proper positioning and avoid lower back stress.

1. Center the breastshield over your nipple.
Note: Do not tilt or overfill containers while pumping. Take care not to kink tubing while pumping.

2. Press the on/off button to start Swing in the Stimulation Phase.

3. Adjust the vacuum to your comfort level while pumping by pressing the increase (+) or decrease (-) button. Holding down these buttons will adjust vacuum more rapidly.
Note: The green LCD light will flash rapidly when you have reached minimum and maximum pumping settings.

4. After two minutes, the Stimulation Phase will automatically change to the Expression Phase. If you let-down (your milk starts to flow) sooner than two minutes, you may pump more efficiently by pressing the let-down button at that time. Adjust to your maximum comfort vacuum™ for most efficient pumping (see next section for details on maximum comfort vacuum).

Maximum Comfort Vacuum™

Research has shown that you will pump more efficiently - meaning get more milk in less time - when pumping at maximum comfort vacuum in the expression phase.

1. Begin in stimulation phase

First, start in the stimulation phase of 2-Phase Expression (this happens automatically in Medela 2-Phase Expression breastpumps, like Swing).

2. Switch to expression phase

Once your milk starts to flow, press the let-down button on the motor to switch to the expression phase. If you don't press the let-down button manually, the expression phase will begin automatically after two minutes.

3. Set your Maximum Comfort Vacuum™

Maximum comfort vacuum is the highest vacuum you can use and still feel comfortable, and is different for every mom. To find your maximum comfort vacuum: increase vacuum until pumping feels slightly uncomfortable (never painful), then decrease slightly.

Mobility

Swing is small and lightweight, and comes with everything you need.

■ Use on any tabletop or flat surface

■ Clip onto your belt or pocket

■ Hang around your neck, over your shoulder, or from a hook

Store your Swing and pump kit (after cleaning and drying) in the included drawstring bag between pumping sessions to keep it clean and all in one place.

Cleaning & Hygiene

Important: *This product is a single user item. Use by more than one user may pose a health risk and performance questions and voids the warranty.*

Hygiene

Wash hands before touching breasts or containers and avoid touching the inside of containers or lids.

Cleaning

Follow these cleaning instructions unless you are told otherwise by your healthcare professional.

Prior to cleaning, disassemble all parts of your pump kit (SoftFit breastshield, valve, membrane, bottles and bottle caps). Be especially sure to remove the yellow valve from the breastshield and separate the white membrane from the valve prior to cleaning. See parts E & F on (page 1).

All parts can be cleaned with soap and water or boiling.

Prior to first use:

Disassemble all parts of your pump kit that come into contact with the breast and milk (SoftFit breastshield, pump body, valve, membrane, bottles and bottle caps). Sanitize by boiling for 10 minutes.

After each use:

Disassemble and wash all parts of your pump kit that come into contact with the breast and milk (SoftFit breastshield, pump body, valve, membrane, bottles and bottle caps) in warm, soapy water.

- Rinse in clear water.
- Air-dry on clean towel and cover parts when not in use.

Do not wash the motor or power cord. Use a damp (not wet) cloth to wipe down the motor if necessary (see part A on page 1)

Alternatively, use Quick Clean™ Breastpump Wipes or Quick Clean™ *Micro-Steam™ Bags for easy, convenient daily cleaning.

Tubing Care:

Inspect tubing after each pumping session for condensation and/or milk.

If condensation appears in tubing:

To eliminate condensation after you have completed pumping, continue running the pump with the tubing attached for another 1-2 minutes or until dry.

*Micro-Steam Bags should not be used to clean Softfit™ Breastshields.

If milk gets into the tubing:

1. Turn off the pump and unplug from power source.
2. Remove tubing from the pump. Wash in soapy water, then rinse in cold clear water.
3. Shake water droplets out and hang to air-dry. For faster drying, attach tubing to the pump and run for 1-2 minutes or until dry, or pour a small amount of isopropyl alcohol through tubing to dry.

If milk accidentally gets inside the motor unit:

1. Turn off the pump and unplug from power source.
2. See troubleshooting instructions, page 11 part C.

Storing Breastmilk

Check with your hospital for specific storage instructions.

When freezing, do not fill containers more than 3/4 full to allow space for expansion. Label containers with the date of pumping.

Breastmilk Storage (For Healthy Term Babies)					
	Room Temperature	Cooler with 3 Frozen Ice Packs	Refrigerator	Self-contained Refrigerator Freezer Unit	Deep Freezer
Freshly expressed breastmilk	4 hours at 66-72°F (19-22°C) ¹	24 hours at 59°F (15°C) ¹	5 - 7 days at 32-39°F (0°C) ²	3-4 months ³	6-12 months at 0°F (-19°C) ³
Thawed breastmilk (previously frozen)	Do not store	Do not store	24 hours ⁴	Never refreeze thawed milk	Never refreeze thawed milk

¹ Hamosh M, Ellis L, Pollock D, Henderson T, and Hamosh P: Pediatrics, vol. 97, No. 4, April 1996, pp 492-497. (4 hours at 77° F/25° C).

² Sosa, Roberto; Barness, Lewis: AJDC, Vol. 141, Jan. 1987.

³ Lawrence R, and Lawrence R: Breastfeeding: A Guide For the Medical Profession, 1999, p.894.

⁴ Arnold L: Recommendations for Collection, Storage and Handling of a Mother's Milk for Her Own Infant in the Hospital Setting, 3rd Edition. The Human Milk Banking Association of North America, Inc. 1999, p. 18.

Important Note: *If the temperature in the room, car or outdoors exceeds 77° F / 25° C, chill milk immediately to preserve freshness.*

Defrosting

- Thaw milk overnight in the refrigerator, or hold the bottle under warm running water to quickly thaw. You can also place the sealed container in a bowl of warm water for 20 minutes to bring it to body temperature.
- Thawed milk is safe in the refrigerator for 24 hours. **Do not refreeze.**

Caution:

- Do not thaw frozen breastmilk in a microwave or in a pan of boiling water.
- Never microwave breastmilk. Microwaving can cause severe burns to baby's mouth from hot spots that develop in the milk during microwaving. Microwaving can also change the composition of breastmilk.
- If adding expressed breastmilk to container of already frozen breastmilk, make sure to add a lesser amount than the already frozen amount.

Troubleshooting

A. If There Is Low or No Suction:

1. Make sure all connections of the kit and pump are secure.
2. Make sure the breastshield forms a complete seal around the breast.
3. Make sure the white membrane is laying flat against the yellow valve. Check that the yellow valve is clean and not damaged. Refer to assembly instructions on page 3.
4. Make sure the tubing end fits securely into the pump body & motor.
5. Take care not to kink tubing. If milk is getting into tubing, check that the tubing is not kinked or pinched.
6. Perform rinsing procedure as described on page 11 (troubleshooting, section C)

B. Pump will not turn on:

1. Make sure pump is plugged in correctly. Pump will NOT turn on if tubing is plugged into the AC adapter port.
2. Make sure the "on" button is being held between 1 and 3 seconds.
3. Do not repeatedly hit the on/off button, or hold the on/off button down longer than 3 seconds.
4. If using batteries, make sure transformer is not plugged into pump.
5. If using batteries, make sure batteries are installed correctly inside pump.

C. If milk accidentally gets into the motor:

Turn off the pump and unplug from the power source. Use a damp (not wet) cloth to wipe the motor and proceed as follows:

1. Remove battery compartment cover from bottom of motor.
2. Remove small, black foam piece from bottom of pump.
3. Rinse the foam piece with clear lukewarm water, and squeeze out excess water. Allow to dry thoroughly before inserting back into pump.
4. Pour about 2 ounces of warm soapy water into a container.
5. Connect the tubing to the motor, leaving the other end unconnected.
6. Connect the power cord to the motor, then plug into electrical outlet (do not power with battery).
7. Hold the end of the tubing that is not connected to the pump in the container with the soapy water.
8. Place the motor over a dry sink and switch on. The water will pump through the tubing into the motor and flow out through the opening beneath the battery cover. Allow it to pump until all the water has been flushed through.
9. Repeat this procedure again with 2 ounces of warm clear water.
10. Leave the motor switched on for an additional 3 minutes so it dries completely.
11. As soon as the foam piece is dry, put it in place again, and re-install the battery compartment cover.
12. If necessary, clean the pump kit according to cleaning instructions.

D. Other troubleshooting:

1. The pump is designed to automatically turn off 30 minutes after last setting adjustment.
2. When traveling outside the U.S., use battery power only. If you experience discomfort while pumping, contact a breastfeeding specialist.

Medela Swing Breastpump

This is a single user product. Use by more than one person may present a health risk and performance questions and voids the warranty.

This product is warranted by Medela, Inc. to the original retail purchaser to be free from defects in material and workmanship for the period of 1 year from the date of purchase. In the event of a defect, Medela will repair or, at Medela's option, replace this product, without charge for such replacement parts or labor. Purchaser shall bear all expense for returning this product to Medela. This warranty does not apply to any product used commercially or which has been subjected to misuse, abuse or alteration.

ANY AND ALL IMPLIED WARRANTIES, INCLUDING THE WARRANTY OF MERCHANTABILITY, ARE LIMITED TO DURATION OF 1 YEAR FROM DATE OF PURCHASE. SOME AREAS DO NOT ALLOW LIMITATIONS ON HOW LONG AN IMPLIED WARRANTY LASTS, SO THE ABOVE LIMITATIONS MAY NOT APPLY TO YOU.

THIS WARRANTY GIVES YOU SPECIFIC LEGAL RIGHTS AND YOU MAY ALSO HAVE OTHER RIGHTS WHICH VARY FROM STATE TO STATE OR IN YOUR COUNTRY.

In the US: Before you make a claim under this warranty, or for additional information, call Medela Customer Service at 800-435-8316, 815-363-1166. You may also contact Medela at customer.service@medela.com.

In other countries: Please contact your local distributor or the location where you purchased your pump; or you may write to Medela at the McHenry, IL, USA address.

**MEDELA, INC. P.O. Box 660
McHENRY, ILLINOIS, 60051-0660 USA
customer.service@medela.com
www.medela.com**

Ordering Parts

If your retailer does not carry the parts you need, you can call 800-435-8316 or 815-363-1166 and order with VISA, MasterCard, American Express, or Discover Card. In other countries than the USA, please contact your local distributor. See our website www.medela.com for a distributor listing.

Use Medela's Breastfeeding National Network (1-800-TELL-YOU) to find a breastfeeding specialist in your area. In the U.S. and Canada, call Medela at 800-435-8316. In other countries, contact your local distributor or the location where you purchased your pump.

Parts (See picture Page 1)

Item	Order Number
a. Motor Unit	#0300011
b. Pump Body	#8007196
c. SoftFit Breastshield	#8100402
d. Tubing	#8007215
e. Valve	#8100475
f. Membrane	#8207013
g. Drawstring Storage Bag	#8000838
h. Shoulder/Neck Strap	#8000734
i. Power Cord	#9207045
j. Bottle	#8117009
k. Bottle Stand	#8100462
l. Bottle Cap	#8107057

Other products sold by Medela	
DayCarrier™ Kit	#57016
Pump & Save™ Bags, 20 Pack	#87033
Pump & Save™ Bags, 50 Pack	#87034
PureLan™ 100, .25 oz. tube	#87021
PureLan™ 100, 1.3 oz. tube	#87022B
Quick Clean Breastpump & Accessory Wipes, 24-pack	#87055
Quick Clean Breastpump & Accessory Wipes, 40-pack individually wrapped	#87059
Quick Clean Micro-Steam Bags	#87024

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